VMDIA Care

The Olympia nursing bed: well-tried, durable and safe. Thanks to several bed panel types, an extensive colour gamut and the possibility to hide the castors from view, the Olympia Care allows you to give a personal touch to the living room.

Carriage

- Very stable high/low adjustment by means of supporting arms.
- Trendelenburg (*) and Reverse Trendelenburg adjustment up to ±17°.
- High-quality design castors Ø 125 mm, with separate brakes.
- High-quality design castors Integral S Ø 125 or 150 mm with central brake (*).
- Design pedal to activate the central brake mechanism (*).
- Hidden central braked castors (Ø 50 mm) are available by means of synthetic cover caps (*).
- The carriage has no actuators, which facilitates its accessibility and its possibility to clean.

Mattress platform

- Three- or four- (*) section mattress platform.
- The backrest is electrically adjustable up to $\pm 70^{\circ}$, with auto-regression.
- Backrest adjustment with intermediate stop at 30° inclination. (if equipped with quick release)
- The knee rest is electrically adjustable up to $\pm 34^{\circ}$ (in case of 4 sections).
- The footrest is manually adjustable to get the Fowler position (in case of 4 sections).
- The legrest is manually adjustable up to $\pm 14^{\circ}$ (in case of 3 sections).
- The electric manipulation of the mattress platform sections is equipped with a clamp safety.
- The mattress platform consists of rounded steel mattress platform profiles.
- The mattress platform sections are equipped with openings for patient restraint straps.
- Telescopic bed extension possible up to 20 cm (not with Pollux sideguards).
- The head end is equipped with 2 cases to install an IV rod or a lifting pole.
- No visible fastening bolts in the foot panel.

Electric bed manipulations

- The electric functions can be activated by means of a hand set or a mini nurse box (*).
- These functions are possible via a key combination on the hand set:
 - Trendelenburg
 - Total lock / unlock (after 3 sec.)

 - personally adjustable transport height
- The hand switch is equipped with a pre-

- programmed sleeping function and button for underbedlight (*)
- A mini nurse box (*) makes it possible to switch off several functions selectively.
- The total lock device (*) blocks all electric functions of the hand switch by means of one single button.
- The optional hand switch (*) can switch off several functions selectively by means of a magnetic key.

Accessories

We would like to refer to our accessory list.

Paneltypes Sideguards

- Midea (standard)
- Valentino (*)
- Puro (*)
- Felina (*)
- Pollux (*)
- Trix (*) Atmosphere (*)
- Twin (*)

(*) Option



Dimensions and technical data - Width: 100,5 cm, Length: 208 cm

- Height variation of the mattress platform:
 - \pm 35,5* cm \pm 78,5 cm (with castors Ø 125 mm)
 - \pm 34,5* cm \pm 77,5 cm (with hidden castors)
- Clearance height of the carriage: L: 170 mm, R: 135 mm (with castors Ø 125 mm) L: 160 mm, R: 125 mm (with hidden castors)
- Safe working load: 250 kg
- Weight (without options): 125 kg
- Matress platform dimensions: 204 x 86 cm
- Mattress dimensions: 195 x 85 cm (*) with Trendelenburg + 2,5 cm



Materials and finishing

- Wood: Beech (massive, veneer), chipboard (melaminated).
- Metal: epoxy-coating RAL 9006 white aluminium.
- Self-lubricating synthetic hinges (= maintenance-free).
- Withstands the common cleaning products.

For more information about the materials and construction: please consult our material lists. We would to refer to our extensive colour gamut.

Options

- 4-section mattress platform.
- Removable solid laminate mattress platform sections.
- Foot extension possible up to + 20 cm.
- Panel types: Puro, Felina or Valentino.
- Bumpers and/or wall spacer.
- Sidequards.
- Total lock device (ACL) or mini-nurse box (ACC).
- Hand switch with selective locking system (incl. Trend).
- Illuminated hand switch and/or nightlight under the bed.
- Control box with Bluetooth.
- Wireless handset for control box with Bluetooth (in case of 4-section mattress platform).
- Wooden surrounds for carriage and/or mattress platform frame
 - Accessory hooks.
- Battery.
- Manual quick release of the backrest (CPR).
- Linen holder.
- Brake buzzer (only with central brake).