

Fero

Wooden relax chair

07578 The Fero relax chair, with lowerable armrests and castors, is equipped with a relax mechanism that has been especially developed for the care sector.

The Classic mechanism can be easily adjusted by pushing yourself off a little from the chair, and this without using any levers.

Support

- 4 wooden laminated legs, connected to each other by means of solid beech support beams.
- Glued mortise and tenon joint.
- The sides are firmly connected to each other by means of a solid wooden cross beam and a metal frame.
- 4 Integral S design castors Ø 125 mm on precision ball bearings, the 2 castors at the back have central brakes and a directional lock.
- Extending, synthetic footrest with a step safety and with an antiskid structure.

Seat

- Removable seat, made of a melaminated MDF board of 12 mm.
- Upholstered with fire-retardant in 2 different density zones to guarantee an excellent seat comfort.
- The rounded seat in front stimulates a fluent blood circulation.
- Comfortable seat height of 53 cm.
- The upholstery is finished with imitation leather and decorative stitching.

Backrest

- Made of an ergonomically preformed plywood shell of 12 mm.
- Upholstered with fire-retardant foam in 3 different density zones to guarantee an excellent seat comfort.
- Equipped with an adjustable and removable head/neck cushion.
- Metal push bar at the back.
- The upholstery is finished with imitation leather and decorative stitching.

Legrest

- Made of a plywood board of 10 mm, upholstered with fire-retardant foam and finished with imitation leather.

Lowerable armrests

- The ergonomic front design of the beech laminated armrests makes it easy to stand up or sit down.
- Standard equipped with dinner tray attachments.
- Two solid chromium-plated axes to adjust the armrests in height in 4 different positions.
- Locked downwards by means of a spring lever.
- In their lowest position, the armrests are on the same level as the seat.

Classic mechanism:

- Easy adjustment in 3 positions, i.e. seat position, reading position, and relax position.
- The relax mechanism first raises the legs, after which it is possible to further incline the entire chair.
- Easy to adjust by just pushing yourself off a little, without needing a lever.
- Guaranteed optimal seat comfort in every position: when the chair is reclined, the angle between the backrest and the seat remains constant.
- Optimal leg support, thanks to the telescopic legrest.
- When the chair is put in the seat position, the legrest is turned inwards so that you can stand up without any problem.

Colours

We would like to refer to our extensive colour gamut.
Sample card + samples available on request.

Accessories


We would like to refer to our
accessorielist.



Dimensions

07578 | Wooden relax Fero:

- L 74 x D 82-146 x H 115 cm
- Seat height: 53 cm
- Seat depth: 51 cm
- Seat width: 53,5 cm
- Height armrests: 73,5-53 cm
- Safe working load: 135 kg



Materials and finishing

- Wood: beech (solid, laminated), Plywood, MDF
- Metal: epoxy coating (light grey), chrome.
- Synthetic.
- Foam: CMHR- and CM-foam (fire-retardant).
- Upholstery: see collection.
- Resistant to the common cleaning products.

For further information about used **materials**, **constructions** and **maintenance**: consult our materialfiles.

Options

- Seat cushion with visco-elastic foam and a bi-elastic cover.
- Swivelling broadened armrests.
- Accessory brackets.
- Restraint bracket.
- Lateral headrest.
- Case for IV rod.
- Backrest with removable headrest.
- Double aluminium design brake pedal.
- Dinner tray support on backrest.

